

# Product Specification Sheet



www.endorphin.net



The Endorphin STSi (Sit to Stand Exerciser) is designed to assist the user in sitting, standing and transferring during hip, knee and total leg physical therapy.

## 330i STS - Sit to Stand Exerciser



Ships via *UPS*  
Shipping Dimensions  
**US / Metric**

85 lbs. / 39 kg  
L - 30" / 76 cm  
W - 22" / 56 cm  
H - 26" / 66 cm

Product Dimensions

70 lbs. / 32 kg  
L - 28" / 71 cm  
W - 46" / 116 cm  
H - 42" / 107 cm

### specifications

**Vendor Part Number:** 330i

**Suggested Retail Price:** \$1895.00

**Function:** The adjustable assistance system has been developed to strengthen muscles, increase range of motion, and provide greater patient stability. The closed-chain assistance load is decreased incrementally until the patient can perform the exercise without assistance.

**Operation:** The therapist no longer has to physically lift and lower the patient as the STS provides the needed support. The assistance is adjusted as a percentage of the patient's body weight and has a maximum user weight capacity of 330 lbs. The therapist is still in control of the patient while the exercise is performed.

**Adjustment:** The seat and support handles are adjustable to the patient's need. The assistance level is selected with a simple adjustment; and adjustable support handles allow the patient to sit and stand in a slow and controlled motion. The seat is designed to stop at 18" from the ground thus conforming to the standard seat height.

**Assistance:** Adjustable assistance allows the patient's progress to be monitored and documented. The support assistance is up to 230 lbs in 5 and 10 lbs increments and is mechanically adjusted to provide a smooth and silent operation.

**Construction:** Durable powder coated 2" 11-gauge tubular steel construction. The STS utilizes state-of-the-art hardware and materials to provide years of uninterrupted service.

**Assembly:** Partial Assembly Required

**Features:** Support hand rails are provided for user stabilization when performing exercise. Wheels are provided to allow the STS to be mobile. This allows the equipment to easily be moved in to a patient's room.



**Endorphin's industry best warranty:**

- 1 year bumper to bumper on all parts
- 5 years on all moving parts
- 20 years on the frame



6901 90th Ave. North, Pinellas Park, FL 33782  
phone 727-545-9848 toll free 800-940-9844 fax 727-546-0613  
www.endorphin.net sales@endorphin.net